# idodontine

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## EN Idodontine acrylic resin. Instructions for use

### Product description

UNIDES A ODI presents IDODENTINE resin, specially indicated for the manufacture of dental pieces, crowns and bridges by means of the muffle curing technique. This is a resin of the IDO series, with a formulation that allows to obtain excellent results working in an easy and practical way. It is available in the A1-D4 and 59-81 colour ranges and stands out for its colour variety, working speed, hardness and colour fastness thanks to its special composition.

### Composition

IDODENTINE powder: PMMA based copolymer > 99% / Catalysts and pigments < 1%.

IDODENTINE liquid: Methyl methacrylate > 90% / Dimethacrylate < 9% / Catalyst < 1%.

### Intended use

IDODENTINE resin allows the fabrication of dental pieces, crowns and bridges.

### Intended users

Dental health professionals, dental technicians and dentists, exclusively.

**Indications** Crowns and bridges.

### Contraindications

Permanent posterior crowns and bridges. Curing of the resin directly in the mouth. Patients with known allergic reaction to any of its chemical components.

### Instructions for use

IDODENTINE is used according to the usual procedure for the production of heat-curing resin crowns and bridges.

**Dosage:** IDODENTINE can be dosed freely at the discretion of the professional, generally the proportion is 12g : 5ml powder/liquid, mixed at room temperature with the aid of a spatula, stirring for 30 seconds.

**Note:** If a more fluid mass is desired, the ratio can be left at 10g : 5ml powder : liquid.

Leave to stand for half a minute and mix again vigorously until a creamy paste is obtained, leave to stand for 15 minutes at room temperature (23°C). The working and curing times depend on the ambient temperature. At 23°C the resin allows potting up to 35 minutes after mixing. Dispense properly in a flask for immersion in cold water. Bring to boil for 30 minutes. Allow to cool slowly.

**Finishing and polishing:** IDODENTINE can be finished and polished by the usual method for PMMA materials.

### Precautions

The product is not resistant to oxidising chemicals and strong acids. During finishing, use a hoover and/or mask to avoid breathing dust. The material is not visible with X-rays. Do not reuse the product and do not use expired materials. Do not vary temperatures and working parameters established by the manufacturer (mechanical modifications of the material could occur, with consequences for the final manufactured product). Do not contaminate the product during the processing phase.

### Interactions

Potential interactions with other products used in the treatment are very unlikely because the product is a fully polymerised acrylic resin-based material. Before using the product, the dentist should consider known interactions and cross-reactions of the product with other materials already present in the patient's oral environment to ensure that the restoration can be placed safely and durably.

### Side effects

In individual cases, rashes, sensitivity or allergies may occur when using the acrylic resin-based prosthesis, in which case do not use it and consult a doctor.

### Risks of reuse

Caution; reuse of used dentures and other dental structures on other persons may lead to infections and transmission of diseases.

### Storage conditions

Store in a cool, dry place away from sunlight. The liquid should be stored at a temperature < 30°C.

### Disposal of the product

The powdered resin and polymerised product do not require specific disposal and can be disposed of in household waste. Monomer residues should be disposed of as chemical waste, do not pour into drains. Contaminated products with biological and/or infection hazards should be disposed of in accordance with applicable regulations.

### Hazard statements / Precautionary statements

**Liquid:** H225 Highly flammable liquid and vapour. H315 Causes skin irritation. H317 May cause an allergic skin reaction. H335 May cause respiratory irritation / P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. P261 Avoid breathing dust/fume/gas/mist/vapours/spray. P262 Do not get in eyes, on skin, or on clothing. P280 Wear protective gloves/protective clothing/eye protection/face protection/hearing protection/... P501 Dispose of container and/or any remaining contents in accordance with applicable regulations. **Powder:** EUH208: Contains (Benzoyl peroxide). May produce an allergic reaction. EUH210: Safety data sheet available on request.

### Information for the patient

The professional should inform the end-user of the following points:

- Normal oral hygiene should be maintained.
  - The material is not visible to X-rays.
  - The material is relatively tough, but strong bites should be avoided.
  - In case of allergic reactions to the material, consult a doctor.
  - Of the following residual risks: Allergies and irritations caused by the material. Breakage and possible ingestion of the product. Deterioration of the product. Incorrect use of the product. Inadequate hygiene and maintenance of the product.
- Manufacturer's responsibility**  
Material for dental use: The manufacturer is not responsible for any damage that may result from use other than in accordance with the instructions for use. Class IIa medical device: Any serious incident involving the device must be reported to the manufacturer and to the competent authority of the Member State in which the incident occurred.

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