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## EN Idopress acrylic resin. Instructions for use

### Product description

UNIDES A ODI presents IDOPRESS resin, especially indicated for the restoration of acrylic removable prostheses, relining and finishing of skeletal prostheses by means of the casting technique. This is a self-curing resin of the IDO series, with a formulation that allows to obtain excellent results working in an easy and practical way. It is available in pink, pink veined and transparent (crystal) colours and stands out for its excellent working performance, hardness and colour stability. IDOPRESS resin is cadmium free.

### Composition

Idopress powder: PMMA based copolymer > 99% / Catalysts and pigments < 1%

Idopress liquid: Methyl methacrylate > 85% / Dimethacrylate < 10% and Catalysts < 5%

### Intended use

Idopress resin enables the repair of removable acrylic prostheses, relining and finishing of skeletons using the casting technique.

### Intended users

Dental health professionals, dental technicians and dentists, exclusively.

**Indications** Removable dental prostheses.

### Contraindications

Fixed dentures in the mouth for an uninterrupted period of more than one month. Patients with known allergic reaction to any of its chemical components.

### Instructions for use

Dosage may be at the discretion of the practitioner. Recommended dosage: **10 g powder and 7 ml liquid** IDOPRESS auto. Mix at room temperature (21-23°C) in a mixing cup and stir for 30 seconds. Its special formulation quickly creates a fluid mass especially suitable for the pouring technique on the model.

### Processing times:

For pouring: 5 minutes at room temperature in pouring conditions.

For packing: The kneading point is reached in 11 - 15 minutes. Before this time it can be modelled or applied with a brush.

**Curing:** In a polymeriser at 2 bar and approx. 50°C for 20 minutes.

**Notes:** Do not store the monomer in places with temperatures below 20°C or in refrigerators, this delays the polymerisation reaction and can also cause the absorption of humidity by the liquid. Thick sprues favour the expansion and the final result of the product.

**Finishing and polishing:** IDOPRESS can be finished and polished by the usual method.

### Precautions

The product is not resistant to oxidising chemicals and strong acids. During finishing use Hoover and/or dust mask to avoid breathing dust. The material is not visible with X-rays. Do not reuse the product and do not use expired materials. Do not vary temperatures and working parameters established by the manufacturer (mechanical modifications of the material could occur, with consequences for the final manufactured product). Do not contaminate the

product during the processing phase.

### Interactions

Potential interactions with other products used in the treatment are highly unlikely because the product is a fully polymerised acrylic resin-based material. Before using the product, the dentist should consider known interactions and cross-reactions of the product with other materials already present in the patient's oral environment to ensure that the restoration can be placed safely and durably.

### Side effects

In individual cases, rashes, sensitivity or allergies may occur when using the acrylic resin-based prosthesis, in which case do not use it and consult a doctor.

### Risks of reuse

Caution; reuse of used dentures and other dental structures on other persons may lead to infections and transmission of diseases.

### Storage conditions

Store in a cool, dry place away from sunlight. The liquid should be stored at a temperature < 30°C.

### Disposal of the product

The powdered resin and polymerised product do not require specific disposal and can be disposed of in household waste. Monomer residues should be disposed of as chemical waste, do not pour into drains. Contaminated products with biological and/or infection hazards should be disposed of in accordance with applicable regulations.

### Hazard statements / Precautionary statements

**Liquid:** H225 Highly flammable liquid and vapour. H315 Causes skin irritation. H317 May cause an allergic skin reaction. H335 May cause respiratory irritation / P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. P261 Avoid breathing dust/fume/gas/mist/vapours/spray. P262 Do not get in eyes, on skin, or on clothing. P280 Wear protective gloves/protective clothing/eye protection/face protection/hearing protection/... P501 Dispose of container and/or any remaining contents in accordance with applicable regulations. **Polvo:** EUH208: Contains (Benzoyl peroxide). May produce an allergic reaction. EUH210: Safety data sheet available on request.

### Information for the patient

The practitioner should inform the end-user of the following instructions:

- Maintain normal oral hygiene and adopt a daily hygiene protocol\*, avoid excessive brushing which may cause abrasive wear, and do not immerse the prosthesis in alkalis, acids or other substances that may damage the acrylic materials.
- At least once a week the denture should be removed and kept for cleaning overnight.
- The material is not visible to X-rays.
- The material is relatively tough, but strong bites should be avoided.
- In case of allergic reactions to the material, consult a doctor.
- Of the following residual risks: Allergies and irritations caused by the material. Breakage and possible ingestion of the product. Deterioration of the product. Incorrect use of the product. Inadequate hygiene and maintenance of the product.

**\*Hygiene protocol:** patients should brush the prosthesis at least once a day; the denture should be removed and gently cleaned. Soak and brush with a soft bristle brush and a non-abrasive denture cleanser to remove food debris, plaque and other deposits.

### Manufacturer's responsibility

Material for dental use: The manufacturer is not responsible for any damage that may result from use other than in accordance with the instructions for use. Class I medical device. Any serious incident involving the device must be reported to the manufacturer and to the competent authority of the Member State in which the incident occurred.

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Unión Dental, S.A.